


INSTRUCTION GUIDE



GPHANTOM LUMBAR

LEARN. PRACTICE. IMPROVE.



INSTRUCTION GUIDE

Thank you for choosing **Gphantom**!

We are a company specialized in the development of solutions for medical training. If you want a specific model for your training, contact us. Here we will help you to handle and optimize the durability of your product.

Your **Gphantom Lumbar** is eligible for our Fidelity Program.

Consult conditions.



LEARN. PRACTICE. IMPROVE.



ABOUT YOUR PRODUCT

The Gphantom Lumbar is a model developed for training procedures in the lumbar spine region. The model has the bone structures from the L2 region up to the sacral hiatus.

With the Gphantom Lumbar it is possible to train ultrasound-guided procedures such as joint injections, infiltrations and anesthesia.


Our needle mark recovery technology allows your Gphantom to recover up to 99% of the needle marks within 48 hours after training, when used with needles with thicknesses equal to or less than 21G. For better use of this needle mark recovery technology, we recommend that you leave your product in a cool place, away from sunlight, for 48 hours after each use; We recommend movements of up to 30 degrees in relation to the product surface, using needles such as the ones mentioned in the previous item. Avoid very long movements (large angles) after needle insertion, as these may shorten the life of your model. In addition, we do not recommend injecting or removing materials from your Gphantom during use.

If you have such a need, please contact us.



ABOUT YOUR PRODUCT

INTERNAL STRUCTURES




B mode ultrasound images of Gphantom Lumbar.

PRODUCT CODE	PBD050514
APPROXIMATE WEIGHT	5,01kg
LENGTH	14cm
HEIGHT	14,5cm
WIDHT	22,5cm



HANDLING AND MAINTENANCE

- ① Always keep your product in a fresh place, sheltered from the sun.



OBS: Expose your Gphantom to temperatures above 30°C (86°F), for long periods, can cause deformations and loss of properties;

- ② For training, remove only the packaging lid, keeping it in the original base;



- ③ After training, clean the product with a paper towel, removing excess gel for ultrasound and, then, wash it under running water without removing it from the base.



Do not store your Gphantom with excess gel, it may cause bad smell, thus reducing product durability.

Do not use detergents or sponges, as these damage the product!



HANDLING AND MAINTENANCE

- ④ After cleaning, store your Gphantom according to instructions contained in item I;
- ⑤ Gphantom Lumbar is a product for educational purposes, intended, **ONLY**, for medical training. In case of doubt, contact us through our relationship channels.

CONTACT US

For any questions or requests, contact us:



contato@gphantom.com.br



+55 1699848-2468